



COOK®

News Release

[News Releases](#)[Photos & Downloads](#)[Contact Public Relations](#)

Important Notice - EAN Bar Code Implementation

With the expansion of global trade and computerization, plain language descriptions of products and services need to be replaced by identification systems that are usable in all trade and industry sectors worldwide.

Beginning April 1, 2003, Cook Incorporated will be implementing the European Article Numbering (EAN) bar code, a uniform and standard numbering system which is currently used throughout the world. This new system will enable the efficient management and distribution of our global products to better serve you.

Of concern and importance is the ability of your system to read this new bar code which is written in a language known as **EAN - Code 128**. Presently, the bar code on Cook products is written in **HIBCC - Code 39**. Please check with your system manager to ensure you will be able to read our new bar code. It would also be appropriate to notify departments who are using current bar code language for inventory management of this change. Our typical alpha-numeric catalog number and Quick Reorder numbers will remain the same on our label at this time.


Manufacturer's Prior HIBCC Bar Code Sticker

COOK®

Cook Incorporated • 750 Daniels Way
P.O. Box 489 • Bloomington IN 47402-0489
Phone (812) 339-2255
U.S. Toll Free (800) 457-4500

PERCUTANEOUS ENTRY NEEDLE

SDN - 18 - 7.0



+ H6980072361E

REORDER#
SDN - 18 - 7.0


*18 GAGE NEEDLE, 7CM LONG
*ACCEPTS UP TO .038" DIA. WIRE GUIDE


INTENDED FOR ONE-TIME USE. STERILE IF PACKAGE IS UNOPENED OR UNDAMAGED. STORE THIS PRODUCT IN A DARK, DRY, COOL PLACE. AVOID EXTENDED EXPOSURE TO LIGHT.

Rx only **STERILE EO**

QUALITY SYSTEM
ISO9001 CERTIFIED

18G X 7CM

LOT NO. **LOT** SAMPLE
USE BY  2008/03
QUICK REORDER# 007236



SAMPLE


Manufacturer's New EAN Bar Code Sticker

COOK®

Cook Incorporated • 750 Daniels Way
P.O. Box 489 • Bloomington IN 47402-0489
Phone (812) 339-2255
U.S. Toll Free (800) 457-4500

PERCUTANEOUS ENTRY NEEDLE

SDN - 18 - 7.0



(01)00827002001664(17)080324(10)SAMPLE

REORDER#
SDN - 18 - 7.0


*18 GAGE NEEDLE, 7CM LONG
*ACCEPTS UP TO .038" DIA. WIRE GUIDE


INTENDED FOR ONE-TIME USE. STERILE IF PACKAGE IS UNOPENED OR UNDAMAGED. STORE THIS PRODUCT IN A DARK, DRY, COOL PLACE. AVOID EXTENDED EXPOSURE TO LIGHT.

Rx only **STERILE EO**

QUALITY SYSTEM
ISO9001 CERTIFIED

18G X 7CM

LOT NO. **LOT** SAMPLE
USE BY  2008/03
QUICK REORDER# 007236



SAMPLE

Manufacturer's New EAN Bar Code System

0 1 0 0 8 2 7 0 0 2 0 0 1 6 6 4 1 7 0 8 0 3 2 4 1 0 S A M P L E
A B C D E F G H

- A. Application Identification - Product Information**
- B. Packaging Indicator**
- C. Product Identification**
- D. Check Digit**
- E. Application Identification - Expiration Date**
- F. Expiration Date (YYMMDD)**
- G. Application Identification - Lot Number**
- H. Lot Number**

[Download New Text Data for Update](#)

[Back to Previous Page](#)

Copyright © Cook Group Incorporated 2003
[Terms of Use](#) | [Privacy Policy](#)